



[www.lillycover.com](http://www.lillycover.com)

For more information on **LILLYCOVER** please scan this QR code with smartphone. The combination of the ambiguous meaning, the flower language of Lily or a woman's name, and the recovery Share the purity and the consistent love with people.

**US Office**  
#3003 N. 1ST ST. #249, SAN JOSE, CA 95134  
[Tel] 408 673 7540  
[Mail] sunhee@lillycover.com

**Korea Office**  
CREATIVE FACTORY 4F, #41260, 471,  
Dongdaegu-ro, Dong-gu, Daegu, Korea  
[Tel] +82 53 756 1140  
[Mail] info@lillycover.com

Search **LILLYCOVER** on SNS!




## How to use **LILLYCOVER** ?



1. Connect the device.



2. Measure the right spot following the description.



3. Press IMP button of LILLYCOVER device.



4. Put IMP sensor over your cheeks.



5. Start self skin diagnosis.



6. Finally, view your result.

Enjoy the analysis caring your skin!

## Application

Empower yourself, your beauty and your confidence from skin diagnosis to skin care solution.

*Do it yourself!*




iOS APP




Android APP




Skin Analysis




Skin History




Plasma Massage



Water Care



UV Care




Moisture Care



**LILLYCOVER**



## Function



### Camera Analysis

- Observe maximum 60X magnification
- Measure oil moisture, pore, flushing, wrinkle, pigment, and sensitive



### Portable and Compact Device

- Take device wherever of whenever it is



### Skin Analysis

- Identify 40 skin types by LILLYCOVER's own algorithm
- Suggest care solution to each type



### Environment Sensing

- Check temperature, humidity, UV index based on the location in real time
- Notify the info of surroundings for skin care



### Plasma Treatment

- Enhance skin absorption
- Promote skin repair and disinfection



### Connected with Application


- Suggest everyday quest to motivate
- Check the history of skin analysis
- Recommend the cosmetics to skin type

## Plasma Massage


Plasma is a material that can penetrate the skin without irritation. You can activate collagen, elastin by putting the device on your skin and massaging. Considering individual skin strength, it can be used in three levels.



## Skin Regeneration




- It helps to increase the absorption rate of cosmetics by 120 times.




- Through plasma that exercise the dermis, it can regenerate cells.

## Relieve Skin



- It can relieve atopic dermatitis and sterilize acne bacteria.



- It helps to improve skin moisturization and relieve dead skin.

Everyday I'm translating your SKIN